



Disposable Hypodermic Needle

Weiyou® Disposable Hypodermic needle is made of high quality stainless steel. Weiyou® Disposable Hypodermic needle is used for injections in various animals. Weiyou® Disposable Hypodermic needle have a low-price advantage and cover most of the European and American markets. We look forward to becoming your long-term partner in China.

Product Description

Weiyou® is a professional manufacturer in China, and we supply Disposable Hypodermic needle which is very high quality. Disposable Hypodermic needle is made of high quality stainless steel. Weiyou® Disposable Hypodermic needle is used for injections in various animals. Weiyou® manufactures and supplies round hub veterinary needles of various sizes, and we has a low-price advantage. Various sizes are available and OEM is acceptable.

Specification

Product name	Disposable Hypodermic needle
Size	14G-24G
Color	Red, yellow, blue, green, etc
Material	Stainless steel, plastic

Features

- Translucent, colour-coded luer lock hub;
- Ultra-sharp, tri-bevelled, anti-coring
- Blister packaging
- Sterile, non-pyrogenic
- Silicone coating assures ease of penetration
- Minimizing animal discomfort





www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com



















Item No.	Needles Gauge and Length		Colour Code	BOX/ Carton
	D×L (G×inch)	D×L (mm×mm)		
NL30716112	16G×1 1/2"	1.6×40	Grey	100//10000
NL307161	16G× 1 "	1.6×25	Grey	100//10000
NL30718112	18G×1 1/2"	1.2×40	Pink	100//10000
NL307181	18G× 1 "	1.2×25	Pink	100//10000
NL3071834	18G×3/4"	1.2×20	Pink	100//10000
NL3071858	18G×5/8"	1.2×15	Pink	100//10000
NL3071812	18G×1/2"	1.2×13	Pink	100//10000
NL30720112	20G×1 1/2"	0.9×40	Yellow	100//10000
NL307201	20G× 1 "	0.9×25	Yellow	100//10000
NL3072034	20×3/4"	0.9×20	Yellow	100//10000
NL3072058	20G×5/8"	0.9×15	Yellow	100//10000
NL3072012	20G×1/2"	0.9×13	Yellow	100//10000
NL30721112	21G×1 1/2"	0.8×40	Green	100//10000
NL307211	21G× 1 "	0.8×25	Green	100//10000
NL3072134	21G×3/4"	0.8×20	Green	100//10000
NL3072158	21G×5/8"	0.8×15	Green	100//10000
NL3072112	21G×1/2"	0.8×13	Green	100//10000
NL30722112	22G×1 1/2"	0.7×40	Bue	100//10000
NL307221	22G× 1 "	0.7×25	Bue	100//10000
NL3072234	22G×3/4"	0.7×20	Bue	100//10000
NL3072258	22G×5/8"	0.7×15	Bue	100//10000
NL3072212	22G×1/2"	0.7×13	Bue	100//10000
NL30723112	23G×1 1/2"	0.6×40	Orange	100//10000
NL307231	23G× 1 "	0.6×25	Orange	100//10000
NL3072334	23G×3/4"	0.6×20	Orange	100//10000
NL3072358	23G×5/8"	0.6×15	Orange	100//10000
NL3072312	23G×1/2"	0.6×13	Orange	100//10000
NL30725112	25G×1 1/2"	0.5×40	Red	100//10000
NL307251	25G× 1 "	0.5×25	Red	100//10000
NL3072534	25G×3/4"	0.5×20	Red	100//10000
NL3072558	25G×5/8"	0.5×15	Red	100//10000
NL3072512	25G×1/2"	0.5×13	Red	100//10000

*Other sizes are also available